

PACS UPGRADE

FEBRUARY 8, 2010

On the morning of Monday, February 8, 2010 from approximately midnight to 6am, the PACS system will be unavailable while being upgraded. The user interface on specialized PACS workstations will change significantly from the “windows-style” folder structure to pre-defined, user selectable work lists (each “folder” will be converted to a work list).

The work lists and user interface make it more efficient to review a series of patient studies and offer more flexibility in defining work lists. For example, a work list for patients from “today” will not reset at midnight, but will include the last 24 hours even when crossing over calendar days.

Staff will be on-site at each facility as the system is made available on February 8. For any questions or assistance during or after the upgrade, please contact the PACS Help Desk at 575-0660

START PAGE OPTIONS

Here are two options to change the appearance of your START PAGE. Neither is required but each enhances functionality. If you make a change you can revert back at any time.

(A) Click PREFS to change the START PAGE to show a list of your patients
(B) Click CONFIGURE to choose filters that contain patients on floors and units at each of the hospitals

Worklist Home Referring Exams Prefs Help

Centricity RIS-IC

Currently signed on: [username]
Last sign-on: 04/21/2009 8:27 AM

Tasks

You have no assigned tasks.

Exams Worklist configure

Click on configure to add items.

CONFIGURE START PAGE – Option A

(A) Click HOME and select REFERRING from drop down list
(B) Click SAVE
(C) Click Close

Do not change any other settings!

User Preferences -- Webpage Dialog

User Preferences help save close

General Passwords Report Creation Desktop Integration

Organization

Default Organization: M-Methodist Hospital

Current Session Organization: M

PACS Viewer

Use Diagnostic Viewer on this Workstation:

Import Imagecast PACS for Radiology 3.5 Personal Filters and Folders: import

Imagecast PACS for Radiology 3.5 Viewer Workstation Location: Main Location

Worklist

Page to automatically launch at Login: Home

Action when the worklist queue is emptied: Go to Worklist

Reset closed exams when queue is emptied:

Related Prior search by: Body Site and Modality

Clinical Info Summary - Initially display: Related Exams

For assistance contact the PACS Helpdesk at 575-0660

CONFIGURE START PAGE – Option B

Select Hospital Floor and Unit Filters

(A) Click on the Filter you want in your list
 (B) Click the arrow to move the filter to the Selected list
 (C) Click the arrow to change the order of the Filters
 (D) Click SAVE

LOCATING PATIENT INFORMATION AND IMAGE VIEW

Options on locating patients

(A) Change to display more patients
 (B) Click on name (blue hyperlink) to bring up exam information
 (C) Click on the yellow arrow or dot to sort any column in ascending/desending order
 (D) Click to find a patient not on your list
 (E) Click to view images

Worklist for [Patient ID] 4 Exams

Additional Search By: (None)

Rows: 1 to 4, of 4 | 10 Rows/Page | Stats on Top

<input type="checkbox"/>	Patient Name	Exam	Sts	Date/Time
<input type="checkbox"/>	[Redacted]	RGEN Bone Density...	C	04/17/2009 9:48 AM
<input type="checkbox"/>	[Redacted]	RGEN Lumbar Spine...	C	04/17/2009 9:30 AM
<input type="checkbox"/>	[Redacted]	RGEN US Breast Bi...	C	04/17/2009 10:57 AM
<input type="checkbox"/>	[Redacted]	RGEN Mammogram Di...	C	04/17/2009 9:55 AM

Clinical Exam Notes

Age at exam: 62 Y, 3 D | Dept No: D00888
 Lumbar Spine - 1 View | Org: BO
 04/17/2009 9:30 AM | Ex. Sts: C | Report Status: Pending Cre | Perf. Resource:

Clinical Info | Memos

This patient has no active allergy records on file.

Diagnostic report text

The report is currently in progress and not available for review.

Signs and Symptoms:

History: H/O BREAST CA, POST MENOPAUSAL
 KM, HOPE, DICE, G, FIRST

Comments: HEALTH, TRICARE, V/O DIAGNOSTIC
 BIL MAMMO, BILATERAL U/S IF
 NEEDED, BONE DENSITY ONE VIEW
 LUMBAR, PREV BO

Pt Loc: BO | Phone: (830) 816-9500

For assistance contact the PACS Helpdesk at 575-0660

IMAGE VIEW

The screenshot shows a web-based image viewer interface. At the top, the title bar reads "Image Viewer -- Webpage Dialog". On the right side of the header, the following information is displayed: "Dt/Tm: 04/17/2009 10:56 AM", "Report Status: Pending Creation", and "S/S:". Below the header is a "Patient History Timeline" section with a dark background. It features a horizontal timeline with several points: "US BREA" (04/17), "MG BREA MUSC" (04/17), "CR SPIN" (04/17), "09:08" (04/17), and "MG BREA" (04/14). A yellow arrow points to the "CR SPIN" entry. Below the timeline, the text "Most Recent: 4/17/2009" and "Least Recent: 4/14/2008" is visible. The main area of the viewer displays a lateral X-ray of a lumbar spine. Above the image, the timestamp "4/17/2009 9:31:45" is shown. Below the image, the text reads "Lumbar Spine - 1 View w/l: 4096/2048", "LATERAL", and "Series #1". A red arrow points from a text box to the close button (X) in the top right corner of the window.

Image Viewer -- Webpage Dialog

Dt/Tm: 04/17/2009 10:56 AM
Report Status: Pending Creation
S/S:

Patient History Timeline

US BREA MG BREA MUSC CR SPIN MG BREA
04/17 04/17 04/17 04/17 09:08 04/14

Most Recent: 4/17/2009 Least Recent: 4/14/2008

4/17/2009 9:31:45

Lumbar Spine - 1 View w/l: 4096/2048
LATERAL
Series #1

Click X to close exam

For assistance contact the PACS Helpdesk at 575-0660